

Estufa - EL Max

Modelos: EL 1.0 - EL 1.1 - EL 1.2 - EL 1.3 EL 1.4 - EL 1.5 - EL 1.6

Modelo	Medidas Internas (alt x Larg x Prof) cm	Medidas Externas (alt x Larg x Prof) cm	Capacidade	Voltagem	Nº de Bandejas	Potência
EL MAX	15,0 x 29,0 x 25,0 cm	30,6 x 40,2 x 30,8 cm	11 lt	Bivolt	02	550 Watts
EL 1.0	23,0 x 25,0 x 23,0 cm	43,0 x 36,2 x 35,5 cm	13 lt	Bivolt	02	550 Watts
EL 2.2	24,0 X 32,0 X 26,0 cm	42,0 x 40,5 x 34,0 cm	22 lt	Bivolt	02	550 Watts
EL 1.1	32,0 X 31,0 X 29,5 cm	52,0 x 42,0 x 39,0 cm	30 lt	Bivolt	03	550 Watts
EL 1.2	34,0 x 34,0 x 34,0 cm	55,0 x 45,0 x 52,0 cm	40 lt	Bivolt	03	550 Watts
EL 1.3	45,0 X 44,5 X 42,5 cm	66,0 x 56,0 x 52,5 cm	85 lt	Bivolt	03	1100 Watts
EL 1.4	53,5 x 48,5 x 42,5 cm	75,0 x 60,0 x 51,0 cm	110 lt	Bivolt	03	1100 Watts
EL 1.5	59,5 x 49,5 x 50,0 cm	84,5 x 61,0 x 58,0 cm	147 lt	Bivolt	03	1100 Watts
EL 1.6	80,0 x 50,0 x 70,0 cm	1,50 x 61,0 x 78,0 cm	280 lt	Bivolt	02	2000 Watts